



# 6000GUI

## MS Windows PC Status Software

### “Resident/Call Management System”

The 6000GUI is a MS Windows based Status-Software that provides facility users with a patient-or room status tracking system. It is designed to work with our Voice-Call products to provide management with PC based Graphical User Interface and statistical data base recording of all transactions for creating printed reports. This system uses our 2590 to interface with the Tech Works network. It provides the user with both “Input” and “Output”. The Graphics are active buttons used either with a Mouse or a “Touch Screen” to give the operator complete control of system functions. Each “Button” can be labeled with the Dial Code and Name or use of the associated point on the system. A total of 512 points can be tracked with this one system. Tabs at the top of the screen indicate multiple networks. Screen tabs can be labeled to indicate a “Wing”, “Floor”, or area of a property or building. This product comes installed on a Tech Works 6000PC to operate properly.

- BENEFITS**
- User Friendly
  - Works on Standard PC
  - Simple to Label
  - Gives current “Status” and Management Reports

### Design Information

Requires  
MS Windows XP-Pro Operating System  
Minimum System requirements:  
Intel Pentium Processor @ 300MHz  
128 Meg RAM  
4 Meg Video RAM  
4 Gig Hard drive space  
PS-2 Mouse Pointer  
1 – RS-232 Serial Port  
VGA 800 x 600 rsi. Minimum  
Monitor not included

## Tech Works™



“Making Specialized Communication Easy”

A Division of: **United Communication Technology, Inc.**  
22349 La Palma Ave., Suite #107, Yorba Linda, CA 92887  
Direct Line: **714-694-1040** Fax: **714-694-1041** Toll Free: **1-800-813-1080**  
Web site **www.tech4people.com**

## **Associated Equipment**

2590 Telephone Interface Master with Caller I.D. and Digital / DTMF Interface / RS232 Universal Network Bridge

## **Video Monitor not included**

## **Architects' and Engineers' Specifications**

A Personal Computer (PC) with a minimum of a 300MHz Processor, 64 Meg RAM, 4 Meg Video RAM, 4 Gig Hard drive space, PS-2 Mouse Pointer, 2 each RS-232 Serial Port, and a VGA Monitor with 800 x 600 rsi. minimum must be provided for the System Status Activity Data Recording and Reporting Software. The PC shall include the Microsoft XP Pro Operating System. The System Status Software shall be a MS Windows based program that provides facility users with a patient-or room status tracking system designed to work with Tech Works Voice-Call products. The program shall provide a Graphical User Interface and statistical database recording of all transactions for creating printed reports. The Graphics shall active buttons used either with a Mouse or a "Touch Screen" to give the operator complete control of system functions. Each "Button" can be labeled with the Dial Code and Name or use of the associated point on the system. A total of 512 points can be tracked with this one system. Tabs at the top of the screen indicate multiple networks. Screen tabs can be labeled to indicate a "Wing", "Floor", or area of a property or building.

The Residence Call PC Software shall be Tech Works Model 6000GUI.